



AIR FLUIDIZED THERAPY BED

The FluidAir™ Air-Fluidized Therapy Bed provides superior pressure reduction for patients that experience pressure ulcers and enable to heal their wounds. The bed is designed to operate in home and institutional settings that have adequate floor structure and support. The FluidAir™ bed provides Air-Fluidized Therapy by supporting the patient on a surface of fluidized microsphere (beads).

A large volume of air is forced upward through the bead mass causing the beads to behave similar to a fluid, helping to reduce interface pressures. The fluidization level of the beads and temperature of the air flowing upward through the bead mass can be adjusted to meet individual patient needs. Fluidization and temperature can be controlled at the main control panel or the patient remote control.



What type of patients benefit from an Air Fluidized Therapy (AFT) bed?

Typically patients with large refractory pressure wounds, posterior skin flaps or burn patients benefit from this therapy. In particular, patients who are unable to routinely turn themselves (paraplegics) will find this therapy beneficial in reducing and healing their pressure wounds. The HydroAire AFT bed is capable of taking on patients weighing up to 300 lbs. and typically patients as tall as 6' 5" in height.

How does the bed work?

The tub of the bed is filled with glass ceramic beads called micro-spheres and the beads actually compose the majority of the weight of the HydroAire AFT bed weighing approximately 900 lbs (bed weight without the beads is 221 lbs). The beads lie on top of a diffuser board. Blowers within the beds pedestals force air across the board and cause the beads to react like a liquid. This act is called fluidization.

What are other advantages of AFT beds?

Because AFT beds have such an effective way of reducing interface pressure, one of the side benefits of the bed is helping patients significantly reduce pain associated with pressure ulcers.

Key Product Features

- Air fluidized (bead-bed) technology.
- Commonly used in the post surgical treatments (skin flaps and skin grafts) .
- Thousands of Micro-sphere (beads) are forced into motion through a high powered blower creating a fluid like surface for pressure re-distribution therapy.
- High air flow/air loss to prevent skin maceration.
- Slant wedge head section, inflates to allow for head elevation for improved respiration.
- Waterproof, vapor permeable top cover.



FEATURES

MAIN CONTROL PANEL

- The main control panel at the foot of the bed contains the main buttons for operating the bed.
- Fluidization Control -Allows user to control the blower speed (full throttle is recommended at most times)
- Temperature Control -Allows the user to control temperature range from 80 to 100 degrees
- Lock-Out Mode- Press this Fluidization Lock-Out mode to disable the fluidization on/off button on the remote hand control, as required
- Note: Fluidization Lock-Out mode should be used at the caregiver's discretion to ensure against unintentional deactivation of fluidization from the Remote Hand Control.

REMOTE HAND CONTROL

- Inflate/Deflate Air Pillows with ON/OFF Button
- On/Off Control - activates or inactivates fluidization

FILTER SHEETS AND COVER SHEET

- Filter sheets (there are two) are the front line barrier between the microsphere beads and the patient. The cover sheet then is placed directly on the filter sheet and is the surface that the patient lies on. Do not pin or clamp items to the filter or cover sheets. Also, cigarette burns, tears caused by sharp objects and pinholes can cause bead leaks. Patient injury, personal injury or equipment damage could occur if beads were to leak through filter and cover sheets.

AIR PILLOWS AND WEDGE CUSHIONS

- The Air Pillows provide upper chest elevation and are controlled with the included Hand Controller. Additionally, Slant Wedge Cushions can be provided on demand. The Wedge Cushion pillowcase can be washed with regular detergent to provide adequate hygiene.



BEADS

- When the FluidAir™ bed is operating the beads become fluid-like and conform to the shape of the patient's body to help relieve pressure. The FluidAir™ bed handles limited amounts of fluids passing through the filter sheet. To contain drainage, use absorbent material or breathable underpads without plastic backing. Excessive incontinence and bodily fluid saturates the beads and hinders fluidization. Petroleum-based topical ointments and silver compounds ruin the coating on the beads and permanently destroy their fluidizing properties. Take care when utilizing these types of products.

SIDE RAILS

- The side rails are intended to be a reminder to the patient of the unit's edges, not a patient-restraining device. Keep the side rail in the up locked position when the patient is in the bed.

DRAINAGE BAG HOLDER

- A drainage bag holder is provided on the side of the FluidAir™ for patients with foley catheters



SAFETY

Unplug Bed - Always unplug the bed immediately after use.

Fluid Spills - Avoid spilling fluids on bed controls. If spills do occur, clean fluids while bed is unplugged while wearing rubber gloves to avoid any possibility of shock. Once fluid is removed, check operation of components affected by the spill.

If this bed has a damaged power cord or plug, is not working properly, has been dropped or damaged, do not operate it. Call the service number located on the bed immediately for service repair.

Power Cord - Keep Power Cord away from heated surfaces.

Smoking - Do not smoke while in the bed.

Air intake - Never block air going into the air intake at the bottom of the blower assembly.

Side Rails - Side rails are designed for the patient's safety. The decision whether and how to use Side Rails lies with the patient and the facility. It is recommended that Side Rails be locked in the full upright position when the patient is unattended. Make sure that the patient knows how to get out of the bed safely (and, if necessary, how to release the Side Rails) in case of fire or other emergency.

Temperature - Recommended room temperature is 80 degrees or lower in a well-ventilated room. Higher room temperatures or poorly ventilated rooms may cause the blower assembly to overheat and shut down, suspending fluidization.

Use - Use this bed only as described in this manual. Do not use attachments not recommended by the manufacturer.

Microsphere Spills - Microspheres (beads) spilled on the floor may be slippery. Any beads on the floor or other exposed surfaces should be wiped up immediately with a damp mop or cloth. After contact with the beads, always wash your hands and avoid rubbing your eyes. If the beads contact the eyes, wash out eyes immediately.

Filter Sheet Condition and Microsphere Containment

- The filter sheets are the essential components in containing beads within the Fluidization Tank. The following guidelines should be observed to ensure filter sheet integrity:

- Do not use filter sheets as draw sheets.
- Do not pin or clip items to filter sheets.
- Objects with sharp edges that must be brought into contact with the filter sheets during patient care should be covered. If an object with sharp edges cannot be covered it should not be used on the surface of the bed.
- Inspect filter sheets regularly for damage, signs of excessive wear or any other evidence of un-serviceability. In situations where the patient is particularly active on the bed surface there is a much greater potential for rips, tears or premature wear.
- Verify that there is sufficient slack in the filter and the Cover sheet to allow patient to be moved freely in any direction.
- Ensure Cover sheet and filter sheets are not wrinkled under patient. Periodically pull both sheets upward to promote better airflow.



GROUNDING INSTRUCTIONS

WARNING: When using this product in the home, use a properly grounded, three-prong, 120 volt outlet. Failure to ground properly could result in personal injury, fire or house wiring.

SHOCK HAZARD: Position the power cord to keep people from tripping over it. When the product is not in use, properly store the power cord away from traffic areas. Failure to do so could result in personal injury.

This product must be grounded. In the event of an electrical short circuit, the grounding reduces the risk of electric shock by providing an escape wire for the electric current. This product is equipped with a grounded power cord and a grounded plug. The plug must be plugged into an outlet that is properly grounded. Close supervision is necessary when this product is used by or near children or invalids.



OPERATING INSTRUCTIONS:

PREPARATION FOR PATIENT PLACEMENT

1. Plug Power Cord into a properly grounded 115 VAC, 15 AMP wall outlet.

NOTE: Do not use a wall outlet controlled by a wall switch.

2. Press the Air On switch on Main Control to turn the FluidAir™ on.
3. Turn Air adjust knob on the Main Control Panel to the Max position.
4. Set Temperature Adjustment on Main Control Panel halfway between Min. and Max. position. Min.=80° Max.= 100°
5. Determine if you wish to Lock-Out the Patient Remote Hand Control by depressing Lock-Out mode in the ON or Off position
6. Utilize a draw sheet over the green Gore-Tex® sheets to assist in positioning the patient.

PATIENT PLACEMENT - PATIENT TRANSFER TO THE FLUIDAIR™ BED

1. Inspect the floor around the bed to verify no beads have spilled.
2. Set the Fluidization level on Main Control Panel to a comfortable position to create a firmer surface.
3. Lower the Side Rails.
4. Place pillow or slant wedge toward the head of the bed.
5. Before placing the patient on the bed, Lay a flat bed sheet on the unit. Do not tie the flat bed sheet down or tuck it in.
6. Transfer patient in accordance with safety rules.
7. Center patient in middle on support surface. Patient's waist should be just below the base of the pillow or slant wedge.
8. Verify there is sufficient slack in the cover and filter sheets to allow patient to be moved freely in any direction. Pulling the sheets too tight can create shear forces for the patient.
9. Ensure that Cover and Filter sheets are not wrinkled under patient. Periodically pull both sheets upward to promote better airflow.
10. Position the Foley Bag Holder as required.
11. Raise the Side Rails.

COMPLETION OF PATIENT PLACEMENT

1. Verify patient is fluidized on the surface properly. Gently rock patient from side to side.
2. Ensure the Cover and Filter sheets are not wrinkled under patient.
3. Press Fluidization Lock-Out mode on the Main Control Panel to disable the Patient Remote Hand Control, as required.
Note: Fluidization lock-Out should be used at the caregiver's discretion to ensure against unintentional deactivation of Fluidization from the Remote Hand Control.
4. Ensure air intake (located at the bottom of the blower) is not blocked by a blanket or other object. The bed has a thermal cutoff that will automatically deactivate the blowers at intervals if the air intake is blocked.