



ROHO® MOSAIC® CUSHION

The Mosaic cushion is especially suited for users who require a cushion for basic level skin protection, support and comfort. Using air flotation technology, the mosaic cushion is not only comfortable, but the use of interconnected **Air** cells to distribute forces evenly across the cushion allows it to be used effectively for low risk skin protection.



Key Product Features

- Lightweight, helps reduce fatigue by adding very little weight to the wheelchair or seating system
- Constructed of polyvinyl, the Mosaic cushion is an affordable solution for skin protection