

AME 3000 LOW AIR LOSS & ALTERNATING PRESSURE MATTRESS



AME 3000 LOW AIR LOSS & ALTERNATING PRESSURE MATTRESS

All air cushions are constructed of high quality and very pliable coated nylon for low shear and even lower interface pressures. Gaymar Plexus Quilted Top Sheet is highly vapor permeable and over-sized providing a low friction, low shear barrier against complications due to moisture.

Machine washable. Unique modular design facilitates cleaning and interchangeability. Convoluted foam base provides additional safety for your patient.



Key Product Features

The Control Unit:

- Effective low-air-loss to assist with moisture control
- Selectable Alternating Pressure with a five minute cycle time
- Max. Inflate for easy patient positioning
- Patient settings are restored after power outages
- CPR quick release provides for rapid deflation
- Control Unit is compact and lightweight

The Therapeutic Support Surface:

- Deep air cells constructed of high quality coated nylon for low shear and effective therapy
- Plexus quilted top sheet is highly vapor permeable and oversized providing a low friction, low shear barrier.
- Machine washable and available in nylon or silver
- Unique modular design facilitates cleaning and interchangeability
- Convoluted foam base provides additional safety for your patient