



AME1500

The American Medical Equipment AME1500 with 3-zoned support surface with vented air cells to provide True Low-Air-Loss Therapy. For added patient protection these air cells rest on a 2" covered foam safety pad. The quilted top cover is waterproof, highly vapor permeable and completely machine washable. The unique CPR release makes deflation quick and easy for correct patient positioning.

The user-friendly control unit is simple to operate with cost effective therapeutic outcomes. The True-Low-Air Loss therapy diffuses under the patient over 100 liters of air per minute. This therapy helps keep the patient dry and comfortable. The adjustable pressure settings are instantaneous for accurate patient placement. The unique filter placement makes hassle free changing or cleaning without using tools or removing filter covers.

Alt-Pressure & Low Air Loss



- Fowler feature for seat inflation when head is elevated
- Alternating cycle time adjusts from 3 to 20 minutes
- Domestic high strength mattress construction for longevity
- Dual modes of therapy
 - Static Low-Air-Loss
 - Alternating Low-Air-Loss
- Lockout to eliminate unintentional settings change
- Low Pressure / Power Failure audible and visual alarms



1-3 alternation



side panel connections



control unit detail

Key Product Features

- Low Air Loss Mattress Replacement System
- 16 nylon-coated, 8" deep air cells provide full immersion to envelope the user
- Low Air loss circulates air for moisture management
- Waterproof, breathable high quality nylon top cover reduces friction and shear
- Max Inflate is used to inflate mattress in seconds
- 2" foam base adds safety and comfort
- CPR quick release for rapid deflation
- Used with raised safety perimeter with egress openings