

## TOTAL LIFT BED™

The Total Lift Bed™ sets a new standard in patients and care-giver treatment safety and comfort. The revolutionary design allows one care giver, using only the pre-programed user friendly handset control, to raise the user from lying to a full standing position and any positions in between with Zero Lift™

The Total Lift Bed™ reduces the caregivers and patient related injuries from getting in and out of bed.

The Concept: **The Total Lift Bed™** is the only bed having a motorized Foot lifter ™, moving in and out, allowing positions which are not possible in any other bed. The Total Lift special Foot Lifter ™ adjusts itself to the user height, and only then will allow the bed to start tilting.





## With just a press of a button you can bring a patient to any position desired.

The Patent-Pending technology of the **Total Lift Bed™** moves the Foot Lifter™ towards the feet of the patient/user, and only after making contact will the bed start the tilting motion, preventing any unwanted movement of the patient in the bed. If positioning the bed to a fully standing position, the Foot Lifter™ will gradually lower the patient so when the bed is in a fully standing position the Foot Lifter™ will be parallel to the floor, allowing the patient to simply walk out of bed, with or without assistance. The patient can also be tilted for treatment, in any desired angle, for as long as needed.











## **Key Product Features**

- Some of the clinical applications the Total Lift Bed:
  - -Hip Replacement.
- -Knee replacement.
- -Spinal Surgery.
- -Cardio Thoracic Surgery.
- -Respiratory problems.
- -Post Abdominal Surgery.
- -Heart and Lungs Transplant.
- Designed to improve safety and quality of life for people who have difficulty performing daily personal care activities, such as getting in and out of bed as well as for the professional nursing staff, caregivers and family members who care for them. Doing More than Any Other Bed: The Vitalgo Total Lift Bed is the ideal solution for safeguarding users and caregivers, such as family members, nurses and therapists.
- The Bed can bring the patient to all possible position, with an easy press on the Handset control.

- Seating with Legs Down.
- Seating with Legs up.
- Treatment position.
- Fully Standing position with foot lifter parallel to floor.
- Unique Smart Break casters design, controlled by the hand set, for maximum stability.

When reaching the appropriate angle, which will avoid body fraction, the Foot Lifter™ will start moving down towards the floor, bringing the patient safely to a standing position with feet on floor level, always keeping his center of gravity to the back. The Total Lift Bed™ excellent for the use of:

-Hospitals

- -Nursing Homes
- -Home care



